



THUASNE



2349 01

ActionReliever
Semi-soft OA knee brace




ActionReliever

Semi-soft OA knee brace



DYNAMIC OFF-LOADING

- Straps creating a dynamic 3-point pressure system: most effective on heel strike
- Combined rigid adjustable hinge on the affected side and flexible stay on the unaffected side
- Self-dosing adjustment possible by tighten the front straps
- For medial or lateral OA




SUSPENSION

- Silicon thread integrated to the knitting
- Cross strap capture system allows the proper anchoring of the straps during the gait cycle

EASE OF USE

- Only 2 steps for fitting



- Effortless magnetic buckles with finger loop: into place automatically
- Anatomical knitting easy to slide up the leg

COMFORT

- Anatomically shaped straps, self oriented and adjustable
- Anatomical knitting with no compression areas at popliteal crease and at brace edges

**INVISIBLE UNDER
MOST CLOTHES**



ActionReliever


Semi-soft OA knee brace

Proven clinically:

Randomised controlled study with 32 patients
(17 bracing group / 15 control group)
on 8 weeks:

Main objective:

Effect of the brace on the pain free walking distance



«By wearing the brace, the pain on exertion was reduced and the pain-free distance extended.»

Benning M, Schneider-Nieskens R. Superiority of a Knee Relief Orthosis in the Treatment of Knee Osteoarthritis. Orthopädie Technik 08/17: 24-30

Medial Left
Lateral Right

Medial Right
Lateral Left

Measurement

Code

Measurement

Code

1	32 - 34 cm	23490120100151	32 - 34 cm	23490120100152	38 - 41 cm	32 - 34 cm	28 - 31 cm
2	35 - 37 cm	23490120200151	35 - 37 cm	23490120200152	42 - 45 cm	35 - 37 cm	32 - 35 cm
3	38 - 41 cm	23490120300151	38 - 41 cm	23490120300152	45 - 48 cm	38 - 41 cm	35 - 38 cm
4	42 - 45 cm	23490120400151	42 - 45 cm	23490120400152	48 - 51 cm	42 - 45 cm	39 - 42 cm
5	46 - 49 cm	23490120500151	46 - 49 cm	23490120500152	50 - 53 cm	46 - 49 cm	41 - 44 cm
6	50 - 53 cm	23490120600151	50 - 53 cm	23490120600152	53 - 56 cm	50 - 53 cm	44 - 47 cm
7	53 - 56 cm	23490120700151	53 - 56 cm	23490120700152	56 - 61 cm	53 - 56 cm	48 - 51 cm
8	56 - 59 cm	23490120800151	56 - 59 cm	23490120800152	61 - 66 cm	56 - 59 cm	51 - 54 cm



THUASNE UK

Unit 4, Orchard Business Centre, North Farm Road, Tunbridge Wells, Kent, TN2 3XF

Tel: 01295 257422, Fax: 01295 25787, Email: orders@thuasne.co.uk

www.thuasne.co.uk