



THUASNE

LombaStab BACK BRACES RANGE

Back braces
with an **easy and precise
tightening⁽¹⁾** that provide
adaptable compression
to the patients.



**MORE EFFICIENCY
AND COMFORT
WITH TECHNICAL
FEATURES**
TO HELP PATIENTS
RETURN BACK
TO THEIR DAILY
ACTIVITIES⁽²⁾



**Automatically adapts
to your body shape** thanks
to the lacing system

More comfort ⁽⁴⁾
thanks to the curved abdominal
panel design and the soft borders
at hips/ribs areas




**Conformable back splint
and height adjustable
to patient's spine**



LombaStab dorso

Height: 35 cm + shoulder straps

Shoulder strap system easy to fit and to adjust



**Strong support
for an efficient
analgesic effect ⁽³⁾**
thanks to the back
plastic plates

Breathability
thanks to the specific elastic
fabric & the Coolmax[®] lining



LombaStab vario

Height: 26 cm

Massaging effect
thanks to the **additional
back pad**



**Easy and precise
tightening ⁽¹⁾**
thanks to the
lacing system



LombaStab

Heights: 21 & 26 cm

LombaStab immo

Height: 35 cm


Removable
and thermoformable
rigid corset



LombaStab high

Heights: 35 & 40 cm


**Easy to use,
ergonomic,
comfortable and
symmetrical
hand loops**
to facilitate the
well-centered
application
of the braces
on the spine





THUASNE

LombaStab range



	Non-specific low back pain (acute, subacute and chronic phases)	Lumbar spinal stenosis	Spondylolysis	Postoperative immobilisation	MODIC type 1 discopathy	Static disorders (kyphosis, lordosis)	Stable traumatic vertebral fracture	Stable osteoporotic vertebral fracture	Scheuermann's disease
LombaStab	⊙	⊙							
LombaStab vario	⊙	⊙							
LombaStab high	⊙	⊙	⊙						
LombaStab immo	⊙	⊙	⊙	⊙	⊙				
LombaStab dorso	⊙					⊙	⊙	⊙	⊙

(1) Study conducted in-house on a panel of 13 people, July 2019

(2) Effects of elastic lumbar belts on the effect of a muscle training program for patients with chronic back pain. Dalichau S, Scheele K. Z Orthop ihre Grenzgeb. 2000 Jan-Feb; 138 (1):8-16. This document is intended solely for the use of healthcare professionals. Lumbar belts are regulated medical devices which bear the CE marking. Read all instructions before use and seek advice from a qualified healthcare professional.

(3) P. Calmels et al. Effectiveness of flexible lumbar support provided by the Thuasne Lombacross Activity support belt in the treatment of subacute low back pain. Spine, 2009; 34(3):215-220 .

(4) Less material at hips and ribs.

Update : January 2020 - Photos credits : ©Thierry Malty - 1905220



THUASNE, 120 rue Marius Auphan 92300 Levallois Perret - France
THUASNE DEUTSCHLAND GmbH - Im Steinkamp 12 · D-30938 Burgwedel